

RED RIVER AUTHORITY OF TEXAS E. coli BACTERIA LOG



Date on:		Time on:			Technician(s):				
Expiration Date of Media:				Start Temp:°C					
	Sample Location	Sample ID No.	mL Use d	Small Cells	Larg e Cells	MPN Table	Dilutio n Factor	MPN/ 100 mL	
1									
2									
3									
4									
5									
6									
7									
8									
9									
10									
11									
12									
13									
14									
15									
16									
17									
18									
19									
20									
Date Off: Ti		Time Off	Time Off:			End Temp:°C			
Method Used to Determine Counts: E. coli IDEXX MPN Chart									
COMMENTS:									
MPN / 100 mL = MPN Table * Dilution Factor QAO									