



# RED RIVER AUTHORITY OF TEXAS

## FECAL COLIFORM BACTERIA LOG



Date on:		Time on:		Technician(s):	
Exp. Date of Media:				Start Temp: °C	
	Sample Location	Sample ID No.	mL Used	Colonies Counted	Fecal Coliform # /100 ml
1	Beginning Blank (#1)	-----			
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20	End Blank (#2)	-----			

Date off:	Time off:	End Temp:	°C	Technician(s):
Method Used to Count Colonies: 20 – 60				
COMMENTS:				
<p>✓ Filter Manifold: ( ) Autoclaved or ( ) Flamed with reagent alcohol prior to use.</p> <p>Fecal Coliform #/100 = (Colonies Counted / mL Used) * 100</p> <p>QAO</p>				